

LABELLING COMPOSITION AND METHOD

5 This invention concerns compositions comprising random mixtures of oligonucleotides and their use for labelling nucleic acids by a random prime method.

Feinberg and Vogelstein (1, 2) introduced the use of random sequence hexanucleotides to prime DNA synthesis on denatured template
10 DNA at numerous sites along its length. The primer-template complex is a substrate for the "Klenow" fragment of DNA polymerase I. By replacing a non-radioactive nucleotide with the radiolabelled equivalent in the reaction mixture, newly synthesised DNA is made radioactive.

Very small amounts of input DNA can be labelled, enabling
15 very high specific activity probes to be produced with relatively small quantities of added nucleotides. These radioactive labelled fragments can then be used as sensitive ^{hybridization} hybridisation probes for a wide range of filter based applications (3-6).

There are several labelling kits that are commercially
20 available for the labelling of DNA by the random prime method. These include the Multiprime, Megaprime, Rediprime and Fluorescein Gene Images kits available from Amersham International plc. Ready-To-Go kits are available from Pharmacia and High Prime kits are available from Boehringer.

25 The Multiprime kit was introduced in the 1980s. It provides different tubes containing the different solutions that enable the user to make up labelling mixtures. One such tube contains a random mixture of 6-mer oligonucleotides, another the polymerase enzyme, and another the supply of nucleotides in the reaction buffer. All these separate solutions
30 are stored frozen at -20°C. The purchaser thaws the different solutions, and adds precise quantities of each to the sample of denatured DNA that is

to be labelled, including a labelled nucleotide. This reaction is then usually incubated at 37°C at which temperature, oligonucleotide annealing and chain extension can occur. However, the reaction may also be incubated at lower temperatures such as an ambient room temperature of about
5 20°C.

The Megaprime kit was introduced commercially in the early 1990s. It is similar to the Multiprime kit, except that 9-mer oligonucleotides are used in place of 6-mers. The Megaprime kit has an advantage over the Multiprime kit, in that 9-mer oligonucleotides anneal more strongly (than do
10 6-mers) to a DNA target and form a hybrid having a higher melting temperature. Thus 9-mers achieve better and more rapid priming of a target than do 6-mers.

The Rediprime kit was introduced commercially in 1994. It comprises a mixture of 9-mer oligonucleotides with a polymerase enzyme and a supply of nucleotides. The mixture is supplied in a freeze-dried
15 state. The freeze-dried mixture also contains a dye for easy visualisation. Dried kits for performing nucleic acid manipulation experiments were described by Ortlepp and McKay in EP 298 669 entitled "Performing nucleic acid reactions". The user reconstitutes the mixture by adding liquid containing the DNA template that is to be labelled, and then liquid
20 containing the labelled nucleotide.

The Ready-To-Go kit was introduced during the 1990s. It is based on a random prime solution containing a random mixture of 9-mer or longer oligonucleotides, which solution is dried by a technique described in
25 EP 383 569. A dye is not present. Like the Rediprime kit, the Ready-To-Go kit can be stored at +4°C or at ambient temperature. Promotional literature emphasises the speed of labelling, which results from the use of 9-mer oligonucleotides.

The High Prime kit is a wet kit containing a random mixture of
30 oligonucleotides. The kit literature does not indicate what length of random oligonucleotides are used, but in the related document EP 649 909 A2, the

- 3 -

use of 6-mer, 9-mer, 12-mer and 15-mer is disclosed. No preferred length of random oligonucleotide is given. The solution is ^{stabilizer} ~~stabilized~~ by the use of glycerol and can be stored at between about -20°C and +4°C.

It can be seen that there has been a trend in commercial kits towards the use of longer oligonucleotides, particularly 9-mers or even longer. Going against this trend, it has been determined by Suganuma, A and Gupta, K C (7) that the use of long random primers, especially 9-mers or longer, reduces the priming efficiency of the random primer reaction. These authors worked on solutions which were used without being dried at any stage. The conclusions of these authors conflict with the findings of the present inventors; which findings are to the effect that, when experiments are done with solutions which are not dried, 9-mers provide more rapid and efficient labelling than do 6-mers, and do not give rise to any problem resulting from self-annealing or self-priming. To the best of applicants' knowledge, the conclusions reported by the authors of (7) have not caused the suppliers of random prime kits to use shorter oligonucleotides.

The present invention is based on the discovery that self-annealing occurs when random 9-mers are used in dried predispensed labelling kits, and that this limits their stability and shelf life. The self-annealing occurs during dispensing and storage when the random 9-mers anneal together to form primer-dimers or primer concatemers. These primer complexes become labelled during the normal labelling reaction, which concomitantly reduces the amount of label that is incorporated into copies of the template that are being synthesised during the reaction. Shorter oligonucleotides are not subject to this problem. The problem is specific to 9-mers (and longer oligonucleotides) used in dried kits.

The invention provides a labelling composition comprising a random mixture of oligonucleotides which are 6-mers to 8-mers, said composition present in a dry state. Preferably the composition also contains at least one of: a polymerase enzyme; a supply of nucleotides for

B chain extension; a labelled nucleotide; a dye; a ^{stabilizer} stabiliser; and a buffer.

As the experimental data below shows, 5-mer oligonucleotides are too short to be useful in dried kits. As the length of the oligonucleotides increases from 6-mers to 9-mers, there is a concomitant
5 increase in the aforementioned self-priming problem. On the other hand, longer oligonucleotides anneal more rapidly and strongly to templates than do shorter ones. Taking into account both these factors, applicants believe that 6-mer oligonucleotides are more preferable than 7-mers which in turn are more preferable than 8-mers.

10 The random mixture of oligonucleotides is present in a dry state. Various drying techniques are possible, including that described in EP 383 569, and also freeze-drying or ^{lyophilization} lyophilisation which is preferred.

B It is possible to use any DNA polymerase enzyme in the labelling reaction, for example Klenow, exonuclease free klenow, DNA
15 polymerase I, T7 DNA polymerase, SequenaseTM, ThermosequenaseTM, so long as the reaction buffer conditions are suitable for the specific enzyme being used.

All four of the nucleotides are preferably present in the composition, whether labelled or unlabelled, and the relative molar
20 concentrations may be adjusted to improve the efficiency of labelling. Also when a labelled nucleotide is present, the equivalent unlabelled nucleotide may also be present to improve the efficiency of labelling, or to control the specific activity of the DNA that is being produced from the labelling reaction.

25 These compositions enable a DNA template to be used to produce copies which are labelled radioactively, for example, by using either phosphate labelled with P-32 or S-35, or by using H-3 or C-14 base labelled nucleotides. Alternatively non-radioactive labels may be used, for example, fluorescein, biotin, digoxigenin, rhodamine and cyanine dyes,
30 may be incorporated when, for example, covalently linked to the base moiety of the nucleotide.

- 5 -

Any ^{stabilizer} stabiliser may be present to protect the activity of the enzyme, for example, trehalose, sucrose, BSA, gelatin. A dye may also be present to allow the dried pellet to be ^{visualized} visualised, before use, and to assist in determining that mixing is thorough.

5 The invention also includes a method of making labelled probes for a nucleic acid template, which method comprises incubating the nucleic acid template under chain extension conditions with the labelling composition as herein described. Preferably the template is DNA. The inventor has found that random 6-mers can give fast labelling kinetics
10 (10 minutes labelling time) by being present at high concentration in the reaction mixture. A preferred concentration is 2-10 O.D./ml in the final reaction with about 5 O.D./ml being most preferable. A probe labelled in this manner is suitable for use in a Southern hybridisation.

All the results shown in the examples show labelling with
15 dCTP-³²P, but this is only as a means to show, and quantitate the amount of self-priming that occurred in each reaction. The reactions are able to label DNA with other labels, both radioactive and non-radioactive, as indicated elsewhere in this specification.

20 References

1. Feinberg, A P and Vogelstein, B, Anal. Biochem., 132: 6-13 (1983).
2. Feinberg, A P and Vogelstein, B, Addendum Anal. Biochem., 137: 266-267 (1984).
- 25 3. Southern, E M, J. Mol. Biol., 98: 503-517 (1975).
4. Thomas, P S, Proc. Nat. Acad. Sci., USA, 77: 5201-5205 (1980).
5. Meinkoth, J and Wahl, G, Anal. Biochem, 138: 267-284 (1984).
- 30 6. Grunstein, M and Hogness, D S, Proc. Natl. Acad. Sci, USA, 72: 3961-3965 (1975).

7. Sugunuma, A and Gupta, K C, Analytical Biochemistry, 224: 605-608 (1995).

Example 1. Manufacture of lyophilised reactions with different random primer lengths:

All primers were diluted to 50 O.D./ml in water. The number of enzyme units was the same in each reaction (7 units).

The amount of each component solution is as follows for a 6 ml scale.

	5 mer reaction mix	6 mer reaction mix	7 mer reaction mix	8 mer reaction mix	9 mer reaction mix
Nucleotide buffer	1.998 ml	1.998 ml	1.998 ml	1.998 ml	1.998 ml
Exo-free Klenow (12 µl) 100 units/µl	1200 units	1200 units	1200 units	1200 units	1200 units
Dilution Buffer	28 µl	28 µl	28 µl	28 µl	28 µl
5 mer primer	1.0 ml				
6 mer primer		1.0 ml			
7 mer primer			1.0 ml		
8 mer primer				1.0 ml	
9 mer primer					1.0 ml
20% Trehalose	1.5 ml	1.5 ml	1.5 ml	1.5 ml	1.5 ml
0.2 mg/ml Xylene Cyanol	0.198 ml	0.198 ml	0.198 ml	0.198 ml	0.198 ml
PF Water	1.264 ml	1.264 ml	1.264 ml	1.264 ml	1.264 ml
Total Volume	6 ml	6 ml	6 ml	6 ml	6 ml

Each reaction mix was dispensed into tubes in 35 µl aliquots, and were freeze dried.

Methods:

1. Nucleotide buffer: Labelling buffer from Nick Translation kit (N5000/N5500 Amersham International plc).
2. Dilution buffer: Storage buffer for enzyme dilution.
- 5 3. Labelling Method: Tubes of DNA for labelling were made up as follows:
5 μ l λ HindIII DNA at 5 ng/ μ l in TE buffer.
40 μ l water.
- 10 Placed all tubes in a boiling water bath (95 to 100°C) for 5 minutes,
placed all tubes on ice for 5 minutes, centrifuged briefly,
then added the denatured DNA solutions to the respective dried reaction tube samples
- 15 added 5 μ l Redivue™ dCTP (α^{32} P) (Product Code AA0005: Amersham International plc) (50 μ l total reaction volume).
Incubated all reactions for 10 minutes at 37°C.
Spotted 2 μ l samples out onto PEI-cellulose tlc plates,
Ran plates in 1.25 M KH_2PO_4 pH 3.4.
- 20 Analysed plates on plate scanner, to measure the %incorporation, %self-priming and %dCTP present at the end of each reaction.
The %self-priming is defined as the % of the total radioactive counts that are situated between the incorporated counts and the counts
25 due to the unincorporated dCTP- ^{32}P .

λ HindIII DNA Labelling with dCTP-³²P (Week 1 Test)

Tube	Primer Type	Tube-1			Tube-2		
		% Incorp	% Self-Prime	% dCTP	% Incorp	% Self-Prime	% dCTP
1, 2	5 mers	62.7	7.9	23.0	54.8	7.7	30.9
3, 4	6 mers	79.9	11.2	2.7	82.1	10.8	2.2
5, 6	7 mers	73.5	17.5	2.8	74.3	15.1	3.7
7, 8	8 mers	68.3	19.1	3.2	65.6	20.4	3.6
9, 10	9 mers	64.9	23.7	3.1	61.5	27.2	3.1

The column headed "% Incorp" shows the percentage of dCTP-³²P incorporated as a chain extension product of a primer- λ Hind III DNA hybrid. The column headed "% Self-Prime" shows the percentage of dCTP-³²P incorporated in a complex involving only primers. The column headed "% dCTP" shows the percent of unincorporated dCTP-³²P. The % dCTP figures were unacceptably high when 5-mer oligonucleotides were used, but were acceptable for 6-mers to 9-mers. Within this range, the % Incorp figures decrease as the oligonucleotide length increases from 6 to 9.

Example 2. Long term stability comparison of dried reactions, nonamers compared with hexamers, 3.5 units of Exo-free Klenow per reaction:

The samples were made up as shown in Example 1, but 6 μ l of Exo-free Klenow was used.

DNA Labelling with dCTP-³²P, results are the averages of the three reactions

Week	Nonamers			Hexamers		
	% Incorp	% Self-Prime	% dCTP	% Incorp	% Self-Prime	% dCTP
3	61.9	17.5	6.2	69.6	9.7	6.3
6	71.4	18.0	4.3	80.8	8.4	4.7
10	65.8	20.2	6.4	75.0	11.9	6.9
16	66.5	16.5	8.0	73.4	11.1	6.5
21	78.3	10.8	3.0	84.3	5.8	2.4
25	42.7	11.6	40.4	55.1	5.4	35.1

5 As these figures show, the % incorporation of dCTP-³²P when using 9-mers was initially lower than when using 6-mers and remained lower on storage of the compositions for up to 25 weeks.

Example 3:

10 Using dried reactions as shown in Example 1, the primer was replaced with water for the reaction drying, and was added later as a separate solution, when the reactions were being used. All reactions were incubated for 10 minutes, and then sampled to measure the % incorporation.

- 10 -

Primer Concentration in reaction O.D./ml	% Incorporation (hexamer primers) Average of two reactions	% Incorporation (nonamer primers) one reaction
6.0	78.3	
5.0	83.2	81.0
4.0	67.7	65.6
2.0	51.5	67.0
1.0	45.1	60.2

It can be seen from these results that the same primer concentration (O.D./ml) is required to achieve the same reaction kinetics, i.e. the same % incorporation in 10 minutes with different random primer lengths. This shows that the molar concentration needs to increase as the primer length is reduced.

Although the above results were obtained using wet reagents, the conclusion would apply also when dry primers are used.

10

Example 4:

Densitometer results of Southern ^{hybridizations} hybridisations

25ng labelling reactions were carried out using the Megaprime Labelling Kit RPN 1606 (Amersham International plc) or using labelled probes from dried nonamer or hexamer labelling reactions made as described above in other examples. Southern blots were ^{hybridized} hybridised for 2 hours at 65°C with the labelled probe under standard conditions and then washed in 2 x SSC, 0.1% SDS, 20 minutes at room temperature, followed by two washes in 0.5 x SSC, 0.1% SDS, for 5 minutes 65°C. The dried blots were detected on X-ray film with 2 intensifying screens and place into a -70°C freezer, for 16 hours. After the film was developed using a film processor it was scanned using a densitometer, then the results were analysed using ImageQuant software.

Kit	Time of test after manufacture	Target	%band intensity of Southern hybridisation cf Megaprime control
9mers	1 week	0.25pg	42.23
9mers	1 week	0.5pg	40.12
9mers	1 week	1.0pg	38.93
6mers	1 week	0.25pg	97.09
6mers	1 week	0.5pg	95.02
6mers	1 week	1.0pg	94.33
6mers	37 weeks	0.25pg	74.58
6mers	37 weeks	0.5pg	80.91
6mers	37 weeks	1.0pg	81.17

Conclusions:

The hexamers used in a dried labelling reaction generate
5 labelled probes which gave a much stronger band intensity than when
nonamers are used, not only when tested initially after 1 week, but even
after an extended period of storage (37 weeks at room temperature).